




Comes with Doggie House Carry Bag



HCPCS
E0570

Available for Private Label
Ask Your Local Drive
Sales Representative



Beagle Pediatric Nebulizer

Item #	Product Description	UOM
18090-BE	Beagle Nebulizer Compressor with Disposable Neb Kit	1/ea
18091-BE	Beagle Nebulizer Compressor with Reusable and Disposable Neb Kit	1/ea

- Child-friendly design to encourage therapy compliance.
- Convenient Doggy House Backpack case comes inside each Beagle Neb box.
- Fun and entertaining doggy stickers included.
- Complete with disposable neb kit, mouthpiece, pediatric aerosol mask, replacement filter and instructions.
- Powerful piston pump is efficient, dependable and quiet.
- Retail packaged.

specifications	
Method of Nebulization	Piston Compressor
Particle Size	0.5 µm to 10 µm
Unit Dimensions	5.6" (W) x 8.6" (H) x 8.6" (D)
Weight	3.1 lbs.
Maximum Pressure	42 PSI
Operating Pressure	13 PSI
Liter Flow	8 lpm
Electrical Requirements	120 VAC 60Hz
Warranty	Limited 5 Year

accessories		UOM
18090F	Replacement Filter for 18070 and 18080 Character Nebs	10/bag
NEB KIT 500	Disposable Neb Kit with Mouthpiece, "T," and 7' Tubing	50/cs
18254	Reusable Neb Kit with Mouthpiece, "T," and 7' Tubing	10/cs
MASK 001P	Pediatric Aerosol Mask	50/cs
22311	Pediatric Dragon Mask	1/ea
18031	Universal Nebulizer Carry Shoulder Bag	1/ea